



America


CERTIFICATE

No. U8 063179 0079 Rev. 00

Holder of Certificate: **STORZ & BICKEL GMBH**

In Grubenäcker 5-9
78532 Tuttlingen
GERMANY

Certification Mark:



Product: **General Medical Devices**

Vaporizer

Model(s): **Mighty + Medic**

Parameters:

Rated input voltage:	5 - 15 VDC
Rated input power:	10- 45 W
Protection class:	internally powered; II
Degree of protection against electric shock:	BF
Degree of water protection:	IP22
Intermittent operation:	7 min. ON / 30 min. OFF

Remarks:
The certificate is valid for USA and Canada.
Other countries might have other requirements,
which were not part of this certification.

Tested according to: ANSI/AAMI ES60601-1:2005/(R)2012
CAN/CSA-C22.2 No. 60601-1:2014

This product was voluntarily tested to the relevant safety requirements referenced on this certificate. It can be marked with the certification mark above. The mark must not be altered in any way. This product certification system operated by TÜV SÜD America Inc. most closely resembles system 3 as defined in ISO/IEC 17067. Certification is based on the TÜV SÜD "Testing and Certification Regulations". TÜV SÜD America Inc. is an OSHA recognized NRTL and a Standards Council of Canada accredited Certification body.

Test report no.: 713181689

Date, 2021-03-26

(Paul Bielawski)