



Ref. Certif. No.

DE 3 - 41162

IEC SYSTEM FOR MUTUAL RECOGNITION OF TEST
CERTIFICATES FOR ELECTRICAL EQUIPMENT (IECEE)
CB SCHEMESYSTEME CEI D'ACCEPTATION MUTUELLE DE
CERTIFICATS D'ESSAIS DES EQUIPEMENTS
ELECTRIQUES (IECEE) METHODE OC**CB TEST CERTIFICATE**
CERTIFICAT D'ESSAI OC

Product

General Medical Devices
Vaporizer System

Name and address of the applicant

STORZ & BICKEL GMBH & CO. KG
In Grubenäcker 5-9
78532 Tuttlingen
GERMANY

Name and address of the manufacturer

STORZ & BICKEL GMBH & CO. KG
In Grubenäcker 5-9, 78532 Tuttlingen, GERMANY

Name and address of the factory

STORZ & BICKEL GMBH & CO. KG
In Grubenäcker 5-9, 78532 Tuttlingen, GERMANY

Ratings and principal characteristics

Rated Input Voltage: 110- 120 V 220-240 V
Rated Frequency: 50-60 Hz
Rated Input Power: 270 VA
Protection Class: II
Applied Part Type: BF
Protection against
Ingress of Water: IP 21
vapormed

Trade mark (if any)

Model/type Ref.

VOLCANO MEDIC 2

A sample of the product was tested and found
to be in conformity withIEC 60601-1:2005
IEC 60601-1:2005/AMD1:2012
IEC 60601-1-6:2010
IEC 60601-1-6:2010/AMD1:2013
IEC 60601-1-11:2015as shown in the Test Report Ref. No.
which forms part of this certificate

713115873

This CB Test Certificate is issued by the National Certification Body
Ce Certificat d'essai OC est établi par l'Organisme **National de Certification**

CB 063179 0073 Rev. 00

Date, 2019-08-23

(Paul Bielawski)



Product Service

Page 1 of 1

TÜV SÜD Product Service GmbH • Certification Body • Ridlerstraße 65 • 80339 Munich • Germany